



Lauden Integrative Pharmacy

1820 41st Ave Suite F

Capitola, CA 95010

Phone (831) 462-9880 Fax (831) 462-9998

Ray M. Reyhani., Compounding Pharmacist

www.LaudenPharmacy.com

DISPENSING INSTRUCTIONS FOR PATIENTS

1 CLICK EQUALS 0.25 GRAM

2 CLICKS EQUALS 0.5 GRAM

[Read all of this information and instructions carefully before dispensing from the Topi-Click™.](#)

Step 1



Pop off the cap of the Topi-Click® container by placing your thumb under the tab on the protective outer cap and pressing upward.

Step 2



Turn the base to dispense the number of clicks specified on label.



When rotating base, watch the Visual Indicators

Indicator Node on Body is fixed at the Click Point



Indicator Column moves with Base. This shows Node in line with column.


LEAVE IN THIS POSITION
AFTER EVERY CLICK
(The Home Position)

Step 3



Apply treatment by rubbing the cream, which is now on the top of the applicator pad, to the area of skin as directed.

Step 4



Prevent evaporation of next dose, by snapping the cap on completely.

Step 5



The **bottom** of the plunger is the indicating line; when it reaches the embossed line **above** the word "Refill," there are approximately 12 clicks left. Approx. 32 clicks remain at the line **below** the word "Refill" as shown above.

The Topi-Click® is a one-time use container and not refillable



PP

Learn more at www.Toni-Click.com